

Application Programming Interface (API)



Table of Contents

Chapter 1	Description	4
1.1. TSI Camera SDK.....		4
1.1.1. Imaging Concepts		4
Chapter 2	C++ API	5
2.1. Enumerations		5
2.1.1. TS_ATTR_ID.....		5
2.1.2. TSI_DATA_TYPE.....		5
2.1.3. TSI_PARAM_FLAGS		5
2.1.4. TSI_CAMERA_STATUS.....		6
2.1.5. TSI_ERROR_CODE		6
2.1.6. TSI_ACQ_STATUS_ID.....		7
2.1.7. TSI_CAMERA_CONTROL_EVENT_ID.....		7
2.1.8. TSI_IMAGE_NOTIFICATION_EVENT_ID.....		7
2.1.9. TSI_PARAM_ID		7
2.1.10. TSI_IMAGE_MODE		9
2.1.11. TSI_HW_TRIG_SOURCE.....		9
2.1.12. TSI_HW_TRIG_POLARITY		9
2.1.13. TSI_OP_MODE.....		9
2.1.14. TSI_EXPOSURE_UNITS.....		10
2.1.15. TSI_ADDRESS_SELECT		10
2.2. Data Types		10
2.2.1. TSI_ROI_BIN		10
2.2.2. TSI_FUNCTION_CAMERA_CONTROL_INFO		10
2.2.3. TSI_FUNCTION_CAMERA_CONTROL_CALLBACK.....		11
2.2.4. TSI_FUNCTION_CAMERA_CONTROL_CALLBACK_EX.....		11
2.2.5. TSI_FUNCTION_IMAGE_NOTIFICATION_CALLBACK.....		11
2.2.6. TSI_FUNCTION_IMAGE_CALLBACK		11
2.3. TsiSdk Class.....		11
2.3.1. Constructor / Destructor		11
2.3.2. Methods.....		12
TsiSdk::Open		12
TsiSdk::Close		12
TsiSdk::GetNumberOfCameras		12
TsiSdk::GetCamera		12
TsiSdk::GetCameraInterfaceTypeStr		13
TsiSdk::GetCameraAddressStr		13
TsiSdk::GetCameraName		13
TsiSdk::ElapsedTime		13
TsiSdk::GetLastErrorStr		13
TsiSdk::GetErrorCode		14
TsiSdk::ClearError		14
TsiSdk::GetErrorString		14
TsiSdk::GetUtilityObject		14
2.4. TsiCamera Class		15
2.4.1. METHODS		15
TsiCamera::Open		15
TsiCamera::Close		15
TsiCamera::Status		15
TsiCamera::GetCameraName		16
TsiCamera::SetCameraName		16

TsiCamera::GetParameter	16
TsiCamera::GetDataTypeSize	16
TsiCamera::SetTextCommand	17
TsiCamera::SetTextCallback	17
TsiCamera::SetParameter	17
TsiCamera:: SetCameraControlCallback	17
TsiCamera:: SetCameraControlCallbackEx	18
The TSI_FUNCTION_CAMERA_CONTROL_CALLBACK_EX type definition above describes the function header for the callback function.	18
TsiCamera:: SetImageNotificationCallback	18
TsiCamera:: SetImageCallback	19
TsiCamera::ResetCamera	19
TsiCamera::FreeImage	19
TsiCamera::FreeAllPendingImages	19
TsiCamera::GetPendingImage	20
TsiCamera::GetLastPendingImage	20
TsiCamera::StartAndWait	20
TsiCamera::Start	20
TsiCamera::Stop	21
TsiCamera::GetAcquisitionStatus	21
TsiCamera::GetExposeCount	21
TsiCamera::GetFrameCount	21
TsiCamera::WaitForImage	22
TsiCamera::ResetExposure	22
TsiCamera::GetErrorCode	22
TsiCamera::ClearError	22
TsiCamera::GetErrorStr	23
TsiCamera::GetLastErrorStr	23
2.4.2. List of Camera Parameters	24
2.5. TsiImage Class	28
2.5.1. DATA MEMBERS	28
2.5.2. METHODS	28
TsiImage::Copy	28
TsiImage::Clone	29
2.6. TsiUtil Class	29
2.6.1. METHODS	29
TsiUtil::WriteImageToTIFF	29
TsiUtil::WriteImageToPNG	29
Chapter 3 Thorlabs Worldwide Contacts	30

Chapter 1 Description

1.1. TSI Camera SDK

The TSI Imaging SDK consists of Libraries, Drivers, Documentation, and Examples that make it easy to incorporate a TSI camera into your system.

1.1.1. Imaging Concepts

Basic steps to get image from camera:

- 1) `TsiSDK *sdk = get_tsi_sdk(0);` //Obtain a reference to the sdk using the static function
- 2) `bool success = sdk->Open();` //Open the SDK
- 3) `int numCameras = sdk->GetNumberOfCameras();` //Get the number of available cameras
- 4) `TsiCamera *cam = sdk->GetCamera(0..numCameras-1);` //Get a reference to a camera
 - a. `cam->Open();` //Open the camera
 - b. Register for desired events (ImageNotification, ControlCallback, etc)
 - c. Set/Get desired camera parameters
 - d. `cam->Start();`
 - e. Acquire images (see step 5)
 - f. `cam->Stop();`
- 5) Access the `TsiImage` from the appropriate callback and use the pointer to the data.
 - a. `TsiImage *image = cam->GetPendingImage();` //Call `GetPendingImage` in Notify Callback
 - i. Alternatively, just receive the image in the parameters of Image Callback
 - ii. Make sure image pointer isn't null
 - iii. Check error codes in callbacks as necessary
 - b. Remember to call `cam->FreeImage(image);` ASAP to clear the buffer for incoming images
- 6) `cam->Close();` //Close the camera when you are done
- 7) `sdk->Close();` //Close the SDK when done
- 8) `release_tsi_sdk(sdk);` //Release the resources used by the SDK

Chapter 2 C++ API

2.1. Enumerations

2.1.1. TS_ATTR_ID

```
typedef enum _TSI_ATTR_ID
{
    TSI_ATTR_NAME,
    TSI_ATTR_DATA_TYPE,
    TSI_ATTR_ARRAY_COUNT,
    TSI_ATTR_FLAGS,
    TSI_ATTR_MIN_VALUE,
    TSI_ATTR_MAX_VALUE,
    TSI_ATTR_DEFAULT_VALUE,

    TSI_MAX_ATTR

} TSI_ATTR_ID, *PTSI_ATTR_ID;
```

2.1.2. TSI_DATA_TYPE

```
typedef enum _TSI_DATA_TYPE
{
    TSI_TYPE_NONE,
    TSI_TYPE_UN8,
    TSI_TYPE_UN16,
    TSI_TYPE_UN32,
    TSI_TYPE_UN64,
    TSI_TYPE_INT8,
    TSI_TYPE_INT16,
    TSI_TYPE_INT32,
    TSI_TYPE_INT64,
    TSI_TYPE_TEXT,
    TSI_TYPE_FP,

    TSI_MAX_TYPES

} TSI_DATA_TYPE;
```

2.1.3. TSI_PARAM_FLAGS

```
typedef enum _TSI_PARAM_FLAGS
{
    TSI_FLAG_READ_ONLY           = 0x00000001,
    TSI_FLAG_WRITE_ONLY          = 0x00000002,
    TSI_FLAG_UNSUPPORTED         = 0x00000004,
    TSI_FLAG_VALUE_CHANGED       = 0x00000008

} TSI_PARAM_FLAGS;
```

2.1.4. TSI_CAMERA_STATUS

```
typedef enum _TSI_CAMERA_STATUS
{
    TSI_STATUS_CLOSED,
    TSI_STATUS_OPEN,
    TSI_STATUS_BUSY,

    TSI_STATUS_MAX

} TSI_CAMERA_STATUS;
```

2.1.5. TSI_ERROR_CODE

```
typedef enum _TSI_ERROR_CODE
{
    TSI_NO_ERROR,
    TSI_ERROR_UNKNOWN,

    TSI_ERROR_UNSUPPORTED,

    TSI_ERROR_PARAMETER_UNSUPPORTED,
    TSI_ERROR_ATTRIBUTE_UNSUPPORTED,

    TSI_ERROR_INVALID_ROI,
    TSI_ERROR_INVALID_BINNING,

    TSI_ERROR_INVALID_PARAMETER_UNDERFLOW,
    TSI_ERROR_INVALID_PARAMETER_OVERFLOW,

    TSI_ERROR_CAMERA_COMM_FAILURE,

    TSI_ERROR_CAMERA_INVALID_DATA,

    TSI_ERROR_NULL_POINTER_SUPPLIED,
    TSI_CAMERA_INVALID_DATA_SIZE_OR_TYPE,

    TSI_MAX_ERROR

} TSI_ERROR_CODE, *PTSI_ERROR_CODE;
```

2.1.6. TSI_ACQ_STATUS_ID

```
typedef enum _TSI_ACQ_STATUS_ID
{
    TSI_ACQ_STATUS_IDLE,
    TSI_ACQ_STATUS_WAITING_FOR_TRIGGER,
    TSI_ACQ_STATUS_EXPOSING,
    TSI_ACQ_STATUS_READING_OUT,
    TSI_ACQ_STATUS_DONE,
    TSI_ACQ_STATUS_ERROR,
    TSI_ACQ_STATUS_TIMEOUT,
    TSI_MAX_ACQ_STATUS_ID
} TSI_ACQ_STATUS_ID;
```

2.1.7. TSI_CAMERA_CONTROL_EVENT_ID

```
typedef enum _TSI_CAMERA_CONTROL_EVENT_ID
{
    TSI_CAMERA_CONTROL_EXPOSURE_START,
    TSI_CAMERA_CONTROL_EXPOSURE_COMPLETE,
    TSI_CAMERA_CONTROL_SEQUENCE_START,
    TSI_CAMERA_CONTROL_SEQUENCE_COMPLETE,
    TSI_CAMERA_CONTROL_READOUT_START,
    TSI_CAMERA_CONTROL_READOUT_COMPLETE,

    TSI_MAX_CAMERA_CONTROL_EVENT_ID
} TSI_CAMERA_CONTROL_EVENT_ID;
```

2.1.8. TSI_IMAGE_NOTIFICATION_EVENT_ID

```
typedef enum _TSI_IMAGE_NOTIFICATION_EVENT_ID
{
    TSI_IMAGE_NOTIFICATION_PENDING_IMAGE,
    TSI_IMAGE_NOTIFICATION_ACQUISITION_ERROR,
    TSI_MAX_IMAGE_NOTIFICATION_EVENT_ID
} TSI_IMAGE_NOTIFICATION_EVENT_ID;
```

2.1.9. TSI_PARAM_ID

```
typedef enum _TSI_PARAM_ID
{
    TSI_PARAM_CMD_ID_ATTR_ID,
    TSI_PARAM_ATTR,
    TSI_PARAM_PROTOCOL,
    TSI_PARAM_FW_VER,
}
```

```
TSI_PARAM_HW_VER ,
TSI_PARAM_HW_MODEL ,
TSI_PARAM_HW_SER_NUM ,
TSI_PARAM_CAMSTATE ,
TSI_PARAM_CAM_EXPOSURE_STATE ,
TSI_PARAM_CAM_TRIGGER_STATE ,
TSI_PARAM_EXPOSURE_UNIT ,
TSI_PARAM_EXPOSURE_TIME ,
TSI_PARAM_ACTUAL_EXPOSURE_TIME ,
TSI_PARAM_HSIZE ,
TSI_PARAM_VSIZE ,
TSI_PARAM_ROI_BIN ,
TSI_PARAM_FRAME_COUNT ,
TSI_PARAM_CURRENT_FRAME ,
TSI_PARAM_OP_MODE ,
TSI_PARAM_CDS_GAIN_INDEX ,
TSI_PARAM_CDS_GAIN = TSI_PARAM_CDS_GAIN_INDEX ,
TSI_PARAM_VGA_GAIN ,
TSI_PARAM_GAIN ,
TSI_PARAM_OPTICAL_BLACK_LEVEL ,
TSI_PARAM_PIXEL_OFFSET ,
TSI_PARAM_READOUT_SPEED_INDEX ,
TSI_PARAM_READOUT_SPEED ,
TSI_PARAM_FRAME_TIME ,
TSI_PARAM_FRAME_RATE ,
TSI_PARAM_COOLING_MODE ,
TSI_PARAM_COOLING_SETPOINT ,
TSI_PARAM_TEMPERATURE ,
TSI_PARAM_QX_OPTION_MODE ,
TSI_PARAM_TURBO_MODE ,
TSI_PARAM_TURBO_CODE_MODE = TSI_PARAM_TURBO_MODE ,
TSI_PARAM_XORIGIN ,
TSI_PARAM_YORIGIN ,
TSI_PARAM_XPIXELS ,
TSI_PARAM_YPIXELS ,
TSI_PARAM_XBIN ,
TSI_PARAM_YBIN ,
TSI_PARAM_IMAGE_ACQUISTION_MODE ,
TSI_PARAM_NAMED_VALUE ,
TSI_PARAM_TAPS_INDEX ,
TSI_PARAM_TAPS_VALUE ,
TSI_PARAM_RESERVED_1 ,
TSI_PARAM_RESERVED_2 ,
TSI_PARAM_RESERVED_3 ,
TSI_PARAM_RESERVED_4 ,
TSI_PARAM_GLOBAL_CAMERA_NAME ,
TSI_PARAM_CDS_GAIN_VALUE ,
TSI_PARAM_PIXEL_SIZE ,
TSI_PARAM_BITS_PER_PIXEL ,
TSI_PARAM_BYTES_PER_PIXEL ,
TSI_PARAM_READOUT_TIME ,
TSI_PARAM_HW_TRIGGER_ACTIVE ,
TSI_PARAM_HW_TRIG_SOURCE ,
TSI_PARAM_HW_TRIG_POLARITY ,
TSI_PARAM_TAP_BALANCE_ENABLE ,
```



```
    TSI_MAX_PARAMS  
  
} TSI_PARAM_ID;
```

2.1.10. TSI_IMAGE_MODE

```
typedef enum _TSI_IMAGE_MODES  
{  
    TSI_IMAGE_MODE_ALLOCATE ,  
    TSI_IMAGE_MODE_STREAM ,  
    TSI_IMAGE_MODE_TRIGGER ,  
  
    TSI_MAX_IMAGE_MODES  
  
} TSI_IMAGE_ACQUISTION_MODES;
```

2.1.11. TSI_HW_TRIG_SOURCE

```
typedef enum _TSI_HW_TRIG_SOURCE  
{  
    TSI_HW_TRIG_OFF ,  
    TSI_HW_TRIG_AUX ,  
    TSI_HW_TRIG_CL ,  
    TSI_HW_TRIG_MAX  
} TSI_HW_TRIG_SOURCE, *PTSI_HW_TRIG_SOURCE;
```

2.1.12. TSI_HW_TRIG_POLARITY

```
typedef enum _TSI_HW_TRIG_POLARITY  
{  
    TSI_HW_TRIG_ACTIVE_HIGH ,  
    TSI_HW_TRIG_ACTIVE_LOW ,  
    TSI_HW_TRIG_POL_MAX  
} TSI_HW_TRIG_POLARITY, *PTSI_HW_TRIG_POLARITY;
```

2.1.13. TSI_OP_MODE

```
typedef enum _TSI_OP_MODE  
{  
    TSI_OP_MODE_NORMAL ,  
    TSI_OP_MODE_PDX ,  
    TSI_OP_MODE_TOE ,  
    TSI_OP_MODE_RESERVED_1 ,  
  
} TSI_OP_MODE, *PTSI_OP_MODE;
```

2.1.14. TSI_EXPOSURE_UNITS

```
typedef enum _TSI_EXPOSURE_UNITS
{
    TSI_EXP_UNIT_MICROSECONDS,
    TSI_EXP_UNIT_MILLISECONDS,
    TSI_EXP_UNIT_MAX

} TSI_EXPOSURE_UNITS, *PTSI_EXPOSURE_UNITS;
```

2.1.15. TSI_ADDRESS_SELECT

```
typedef enum _TSI_ADDRESS_SELECT
{
    TSI_ADDRESS_SELECT_IP,
    TSI_ADDRESS_SELECT_MAC,
    TSI_ADDRESS_SELECT_ADAPTER_ID,
    TSI_ADDRESS_SELECT_MAX

} TSI_ADDRESS_SELECT;
```

2.2. Data Types

2.2.1. TSI_ROI_BIN

```
typedef struct _TSI_ROI_BIN {

    UNS32 XOrigin; // X Origin in sensor array
    UNS32 YOrigin; // Y Origin in sensor array
    UNS32 XPixels; // Frame width in pixels
    UNS32 YPixels; // Frame height in pixels

    UNS32 XBin;    // Binning factor in the X dimension
    UNS32 YBin;    // Binning factor in the Y dimension

} TSI_ROI_BIN, *PTSI_ROI_BIN;
```

2.2.2. TSI_FUNCTION_CAMERA_CONTROL_INFO

```
typedef struct _TSI_FUNCTION_CAMERA_CONTROL_INFO {

    uint32_t FrameNumber;

    struct {

        uint32_t Year;
        uint32_t Month;
        uint32_t Day;
        uint32_t Hour;
        uint32_t Min;

    }
```

```
uint32_t Sec;  
uint32_t MS;  
uint32_t US;  
  
} TimeStamp;  
  
} TSI_FUNCTION_CAMERA_CONTROL_INFO,  
*PTSI_FUNCTION_CAMERA_CONTROL_INFO;
```

2.2.3. TSI_FUNCTION_CAMERA_CONTROL_CALLBACK

```
typedef void (*TSI_FUNCTION_CAMERA_CONTROL_CALLBACK)  
(int ctl_event, void *context);
```

2.2.4. TSI_FUNCTION_CAMERA_CONTROL_CALLBACK_EX

```
typedef void (*TSI_FUNCTION_CAMERA_CONTROL_CALLBACK_EX )  
(int ctl_event, TSI_FUNCTION_CAMERA_CONTROL_INFO  
*ctl_event_info, void *context);
```

2.2.5. TSI_FUNCTION_IMAGE_NOTIFICATION_CALLBACK

```
typedef void (*TSI_FUNCTION_IMAGE_NOTIFICATION_CALLBACK)  
(int notification, void *context);
```

2.2.6. TSI_FUNCTION_IMAGE_CALLBACK

```
typedef void (*TSI_FUNCTION_IMAGE_CALLBACK)  
(TsiImage *tsi_image, void *context);
```

2.3. TsiSdk Class

This class provides a framework for all the other classes in the API to exist within.

Instantiating an object of this type and then calling the OpenSDK method instructs the API to examine your system, searching for and building a list of cameras that it finds.

The CloseSDK (coupled with a subsequent call to OpenSDK) can be used to handle instances where the list of cameras changes dynamically.

The destructor for the TsiSdk class will first close and delete any camera objects that were allocated, and then close and clean up the SDK.

2.3.1. Constructor / Destructor

The constructor for the TsiSDK class initializes the library but does not actually build a list of cameras – this functionality is implemented in the OpenSDK method so that the TsiSdk class does not have to be deleted and recreated when the list of cameras changes.

The destructor for this class contains an implicit call to the CloseSDK (see below) method.

2.3.2. Methods

TsiSdk::Open

```
bool Open();
```

This method must be called before any other

TsiSdk::Close

```
bool Close();
```

This method iterates over the list of cameras, first closing any cameras that have been left open, then deleting their objects. After calling this method (or destructing an object of type TsiSdk), make sure your program does not access any TsiCamera objects.

TsiSdk::GetNumberOfCameras

```
int GetNumberOfCameras ();
```

This method is used after the OpenSDK method to get a count of the number of cameras the OpenSDK method found.

TsiSdk::GetCamera

```
TsiCamera* GetCamera (int CameraNumber);
```

This method fails if the value returned is NULL.

Otherwise, the pointer returned is to a camera object of type TsiCamera. See TsiCamera class definition below for more details.

TsiSdk::GetCameraInterfaceTypeStr

```
char* GetCameraInterfaceTypeStr (int camera_number );
```

This method returns a string describing the interface type of the camera.

TsiSdk::GetCameraAddressStr

```
char* GetCameraAddressStr (int camera_number,  
                           TSI_ADDRESS_SELECT address_select );
```

This method returns the specified address type for a specified camera.

TsiSdk::GetCameraName

```
char* GetCameraName (int camera_number );
```

This method returns a string with the name of the selected camera.

TsiSdk::ElapsedTime

```
uint64 ElapsedTime (uint64 StartTime);
```

This method provides a convenient way to calculate the elapsed time in milliseconds between successive calls to this method.

Calling this method with a start time of zero establishes the start time. Calling the method again, using the results of the first call will result in the elapsed time in milliseconds between those calls being returned.

TsiSdk::GetLastErrorStr

```
char* GetLastErrorStr();
```

Returns the string describing the last error.

TsiSdk::GetErrorCode

```
TSI_ERROR_CODE GetErrorCode();
```

Returns the current error code for the SDK, a returned value of zero indicates success. Clears any non-zero error code.

TsiSdk::ClearError

```
bool ClearError ();
```

Clears any existing error code for the SDK without having to read it using GetErrorCode.

TsiSdk::GetErrorString

```
bool GetErrorString(TSI_ERROR_CODE ErrorCode, char  
*StringBuffer, int &StringLength)
```

Returns the non-localized ASCII string description of the supplied error code.

Passing a NULL pointer and pointer to a variable containing a length of zero length results in the length of the string associated with the supplied error code being returned so that a buffer of appropriate length can be dynamically allocated.

TsiSdk::GetUtilityObject

```
TsiUtil* GetUtilityObject ();
```

Returns a pointer to an instance of TsiUtil which has some functions for saving out 16-bit tif and png image formats.

2.4. TsiCamera Class

Objects of this class are used to control and acquire data from the cameras in your system. You cannot simply instantiate an object of this type. Instead, instantiate an object of the **TsiSdk** class, call the **OpenSDK** method, and then use the **GetCamera** method in that class to get a pointer to an object of this type.

2.4.1. METHODS

TsiCamera::Open

```
bool Open ();
```

Opens this camera for subsequent communications. The reason for opening and closing cameras is to conserve resources for cameras that have been enumerated, but are not currently in use.

This method must be called (once) before any other methods other than GetCameraName can be used.

TsiCamera::Close

```
bool Close ();
```

This method is called implicitly by the TsiSdk method CloseSDK and its destructor, but one can call this method directly to close the camera to conserve resources when a camera is not needed.

TsiCamera::Status

```
bool Status (TSI_CAMERA_STATUS *CameraStatus);
```

This method returns the current status of the camera (see the definition of TSI_CAMERA_STATUS above). The intent of this method is to return the overall state of the camera (open, closed, busy, etc.).

For the status of an acquisition, use the GetStatus method of the TsiAcquisition object below.

TsiCamera::GetCameraName

```
char* GetCameraName ();
```

This function returns a pointer to the camera name. If no name has been previously assigned to the camera, a generic camera name will be returned.

TsiCamera::SetCameraName

```
bool SetCameraName (char *CameraName);
```

This function allows the caller to give the camera a user generated name. This function is intended for use in multi-camera systems where the user may want to assign a name to a particular camera that may have a particular function or be connected to a particular part of their system. This name is persistent and is tied to the model and serial number of the camera.

TsiCamera::GetParameter

```
bool GetParameter (TSI_PARAM_ID ParameterID, size_t Length, void *data);
```

This method is the main way of querying the various camera parameters, discovering whether a particular setting is supported on a camera, and other attributes about the parameter, such as its data type, length (in the case of vector or array type), and any other flags the parameter might have.

To use this method, you must supply a parameter identifier (see `TSI_PARAMETER_ID` as defined above), the size of the variable to hold the data, and a pointer to the variable to hold the data returned.

Attributes **about** the parameter, with the exception of the minimum and maximum values a parameter may have all have fixed sizes. The actual value, the min and the max values will all have the same type, and that type can be discovered by using the **TSI_ATTR_DATA_TYPE** attribute identifier. The value returned will be a **TSI_DATA_TYPE** enumerated type. The `GetDataTypeSize` method can be used to determine the scalar size for the type. The total storage required for a parameter will be the size for the data type multiplied by the value returned by **TSI_ATTR_DATA_COUNT** attribute id.

TsiCamera::GetDataTypeSize


```
int32 GetDataTypeSize (TSI_DATA_TYPE DataType);
```

This method returns the size in bytes of any of the SDK defined data types.

TsiCamera::SetTextCommand

```
bool SetTextCommand (char *str);
```

This method sends a serial text command directly to the camera.

TsiCamera::SetTextCallback

```
bool SetTextCommand (TSI_TEXT_CALLBACK_FUNCTION func, void  
*context);
```

Using this method, you can register for a callback to receive serial text responses to serial text commands issued via SetTextCommand(char *str).

TsiCamera::SetParameter

```
bool SetParameter (TSI_PARAM_ID ParameterID, void *data);
```

This method is a sister function to the GetParameter function described above and is used to set any parameter that is does not have the **TSI_FLAG_READ_ONLY** flag set.

A parameter's setting may be cached (unless the API determines it's value is volatile and must be sent to the camera each time, even if the value does not actually change). This is used to cut down on unnecessary traffic to / from the camera.

TsiCamera:: SetCameraControlCallback

```
bool SetCameraControlCallback (
```

```
TSI_FUNCTION_CAMERA_CONTROL_CALLBACK    func ,  
void *                                  context  
);
```

This method allows the caller to specify a callback function for interesting events about the camera.

See the **TSI_CAMERA_CONTROL_EVENT_ID** enumeration above for the supported list of events.

The **TSI_FUNCTION_CAMERA_CONTROL_CALLBACK** type definition above describes the function header for the callback function.

TsiCamera:: SetCameraControlCallbackEx

```
bool SetCameraControlCallbackEx (  
    TSI_FUNCTION_CAMERA_CONTROL_CALLBACK_EX    func ,  
    void *                                      context  
);
```

This method allows the caller to specify a callback function for interesting events about the camera.

See the **TSI_CAMERA_CONTROL_EVENT_ID** enumeration above for the supported list of events.

The *TSI_FUNCTION_CAMERA_CONTROL_CALLBACK_EX* type definition above describes the function header for the callback function.

TsiCamera:: SetImageNotificationCallback

```
bool SetImageNotificationCallback (  
    TSI_FUNCTION_IMAGE_NOTIFICATION_CALLBACK    func ,  
    Void *                                      context  
);
```

Using this method, you can register callbacks for events such as cameras being unplugged or other asynchronous events occurring on the camera.

See the **TSI_IMAGE_NOTIFICATION_EVENT_ID** enumeration above for the supported list of events.

The **TSI_FUNCTION_IMAGE_NOTIFICATION_CALLBACK** type definition above describes the function header for the callback function.

TsiCamera::SetImageCallback

```
bool SetImageCallback (  
    TSI_FUNCTION_IMAGE_CALLBACK      func,  
    Void *                           context  
);
```

This method allows the caller to specify a callback function that will be called whenever there is a new frame of image data available.

The **TSI_FUNCTION_IMAGE_CALLBACK** type definition above describes the function header for the callback function.

TsiCamera::ResetCamera

```
bool ResetCamera ();
```

This method is used to set the camera to its power-on state. The application will have to query the camera parameters to determine what the current camera settings are.

This can be accomplished by iterating over all the camera settings calling a `GetParameter` with the Parameter IDs the application cares about, specifying the `TSI_ATTR_FLAGS` attribute, and checking for the `ATTR_VALUE_CHANGED` flag. Any parameter that has changed as a result of restoring the camera to its power on state will have this flag set.

TsiCamera::FreeImage

```
bool FreeImage (TsiImage *Image);
```

This method frees the image data object. The `Close` method and the destructor for this class will close and free up any image data objects that have been created. Any pointers to `TsiImage` objects that have been created must not be accessed after the object has been freed by the `TsiCamera FreeImageData` or `Close` methods or the `TsiSdk` object's `CloseSDK` method or destructor.

TsiCamera::FreeAllPendingImages

```
bool FreeAllPendingImages (TsiImage *Image);
```

This method frees all of the the pending images in the buffer.

TsiCamera::GetPendingImage

```
TsiImage* GetPendingImage (void);
```

This method will return a pointer to the next available image.

TsiCamera::GetLastPendingImage

```
TsiImage* GetLastPendingImage (void);
```

This method will return a pointer to the last available image.

TsiCamera::StartAndWait

```
bool StartAndWait (int32 Timeout);
```

This method starts an aquisition of one or more frames using the current camera parameters and waits for it to finish. If the elapsed time to acquire the image data exceeds the timeout value (in milliseconds), the acquisition will be cancelled and the function will return FALSE.

Also, if the combination of current camera parameters do not permit proper operation of the camera, FALSE will be returned. An example would be calling StartAcquisitionAndWait with the TSI_PARAM_FRAME_COUNT parameter set to 0 (which indicates acquire until a StopCamera call is made).

This method is intended for use in simple applications that could be simple command line applications that are going to be invoked periodically.

TsiCamera::Start

```
bool Start ();
```

This method starts an acquisition, but does not block.

The application may specify a callback function that will be called as each frame is completed and ready for retrieval using the ReadXxxImage functions.

The application must call the StopCamera method before issuing another StartAcquisition call.

TsiCamera::Stop

```
bool Stop ();
```

This method is used to either cancel and / or clean up after an acquisition. This method must be called between subsequent calls to the StartCamera method.

TsiCamera::GetAcquisitionStatus

```
TSI_ACQ_STATUS_ID GetAcquisitionStatus ();
```

Returns the current status of the acquisition for applications that wish to poll for completion.

TsiCamera::GetExposeCount

```
int32 GetExposeCount ();
```

The GetExposeCount method returns the number of exposures that have been completed since a StartCamera call.

TsiCamera::GetFrameCount

```
int32 GetFrameCount ();
```

The `GetFrameCount` method returns the number of frames that have been completely transferred to host memory since a `StartCamera` call.

TsiCamera::WaitForImage

```
bool WaitForImage ();
```

An application that wishes to poll for frame completion can call this function that will return `TRUE` when the next frame is available.

TsiCamera::ResetExposure

```
bool ResetExposure ();
```

If the camera is waiting for a trigger or in the middle of a long exposure, this method can be used to dump the charge on the sensor, and reset the camera state to the beginning of the exposure.

TsiCamera::GetErrorCode

```
TSI_ERROR_CODE GetErrorCode();
```

Returns the current error code for the camera, a returned value of zero indicates success. Clears any non-zero error code.

TsiCamera::ClearError

```
VOID ClearError ();
```

Clears any existing error code for the camera without having to read it using `GetErrorCode`.

TsiCamera::GetErrorStr

```
bool GetErrorStr(TSI_ERROR_CODE ErrorCode, char *StringBuffer,  
int32 &StringLength)
```

Returns the non-localized ASCII string description of the supplied error code.

Passing a NULL pointer and pointer to a variable containing a length of zero length results in the length of the string associated with the supplied error code being returned so that a buffer of appropriate length can be dynamically allocated.

TsiCamera::GetLastErrorStr

```
char* GetLastErrorStr (void)
```

Returns a pointer to the non-localized ASCII string description of the current error code.

2.4.2. List of Camera Parameters

This list of parameters are used by the **TsiCamera.GetParameter** and **TsiCamera.SetParameter** methods.

TSI_PARAM_ATTR	TSI_TYPE_UNUS8: For expert use only, contact us for more info.
TSI_PARAM_PROTOCOL	TSI_TYPE_UNUS32: Obsolete
TSI_PARAM_FW_VER	TSI_TYPE_TEXT: Returns the firmware revision of the camera. The format of the returned data is MM.mm.rr.build, where MM is a two digit major version, mm is a two digit minor version, rr is a two digit revision, and build is a sequential number used internally by TSI supplied utilities.
TSI_PARAM_HW_VER	TSI_TYPE_TEXT: Returns the hardware rev of the camera. The format of the returned data is MM.mm.rr.build, where MM is a two digit major version, mm is a two digit minor version, rr is a two digit revision, and build is a sequential number used internally by TSI supplied utilities.
TSI_PARAM_HW_MODEL	TSI_TYPE_TEXT: Returns the hardware model name of the camera, format TBD.
TSI_PARAM_HW_SER_NUM	TSI_TYPE_TEXT: Returns the serial number of the camera, format TBD.
TSI_PARAM_CAMSTATE	TSI_TYPE_UNUS32 (enum): Get the current state of the camera (IDLE / RUN)
TSI_PARAM_CAM_EXPOSURE_STATE	
TSI_PARAM_CAM_TRIGGER_STATE	TSI_TYPE_UNUS32 (enum): Get the camera's current trigger state (NONE/ARMED/DISARMED/TRIGGERED)
TSI_PARAM_EXPOSURE_UNIT	TSI_TYPE_UNUS32 (enum): Set the camera's exposure time unit of measurement.
TSI_PARAM_EXPOSURE_TIME	TSI_TYPE_UNUS32: Set the camera's exposure time in the current unit of measurement. See TSI_PARAM_EXPOSURE_UNIT above.
TSI_PARAM_ACTUAL_EXPOSURE_TIME	TSI_TYPE_UNUS32: Gets the actual exposure time in microseconds for a commanded exposure time.
TSI_PARAM_FRAME_TIME	TSI_TYPE_UNUS32: Get the camera frame readout time in milliseconds.
TSI_PARAM_VSIZE	TSI_TYPE_UNUS32: Get the camera sensor's number of active rows.
TSI_PARAM_HSIZE	TSI_TYPE_UNUS32: Get the camera sensor's maximum number of active columns.

TSI_PARAM_ROI_BIN	TSI_TYPE_UN32 [6]: Get/Set ROI and Binning Factors – ROI and binning are returned as an array of six UNS32 values, which can be encoded/decoded using the \$\$\$ structure type defined above.
TSI_PARAM_FRAME_COUNT	TSI_TYPE_UN32: Sets number of frames to acquire before going to the TSI_ACQ_DONE state.
TSI_PARAM_CURRENT_FRAME	TSI_TYPE_UN32: Gets the number of frames completely read out in an "N" frame sequence.
TSI_PARAM_OP_MODE	TSI_TYPE_UN32 (enum): Sets the current operating mode
TSI_PARAM_CDS_GAIN	TSI_TYPE_UN32 (enum): Gets or Sets CDS gain setting
TSI_PARAM_VGA_GAIN	TSI_TYPE_UN32: Gets or Sets VGA gain
TSI_PARAM_GAIN	TSI_TYPE_UN32: Gets or Sets linear gain (x1 – xN)
TSI_PARAM_OPTICAL_BLACK_LEVEL	TSI_TYPE_UN32: Sets the optical black level
TSI_PARAM_PIXEL_OFFSET	TSI_TYPE_UN32: The offset for each pixel
TSI_PARAM_READOUT_SPEED_INDEX	TSI_TYPE_UN32: Selects the index into the set of readout speeds the camera supports.
TSI_PARAM_READOUT_SPEED	TSI_TYPE_UN32: Returns the digitization rate in MHz for the given speed index.
TSI_PARAM_FRAME_RATE	TSI_TYPE_UN32: Returns the frame rate in frames per second. Value is returned as one 32 bit number. The upper 16 bits represent the whole number of frames per second. The lower 16 bits represent the fractional part in hundredths of frames per second.
TSI_PARAM_COOLING_MODE	TSI_TYPE_UN32 (enum): Sets the cooling mode.
TSI_PARAM_COOLING_SETPOINT	TSI_TYPE_INT32: Sets the current CCD chamber temperature set point in hundredths of degrees centigrade (if supported by the camera). This value contains an implied decimal point, to convert to floating point degrees centegrade, divide by 100.0. Inversely, to convert from floating point degrees celcius, multiply that value by 100.0 and round or truncate (cast) as desired.
TSI_PARAM_TEMPERATURE	TSI_TYPE_INT32: Gets the current CCD chamber temperature (if supported by camera). Uses the same format / data type as TSI_PARAM_COOLING_SETPOINT.
TSI_PARAM_QX_OPTION_MODE	TSI_TYPE_UN32 (enum): Sets the QX mode
TSI_PARAM_TURBO_CODE_MODE	TSI_TYPE_UN32 (bool): Sets TURBO Code Output On/Off
TSI_PARAM_XORIGIN	TSI_TYPE_UN32: Gets or sets the X component of the upperleft corner of the specified ROI.

TSI_PARAM_YORIGIN	TSI_TYPE_UN32: Gets or sets the Y component of the upperleft corner of the specified ROI.
TSI_PARAM_XPIXELS	TSI_TYPE_UN32: Gets or sets the width in pixels of the desired ROI.
TSI_PARAM_YPIXELS	TSI_TYPE_UN32: Gets or sets the height in pixels of the desired ROI.
TSI_PARAM_XBIN	TSI_TYPE_UN32: Used to set or query the current horizontal binning factor.
TSI_PARAM_YBIN	TSI_TYPE_UN32: Used to set or query the current horizontal binning factor.
TSI_PARAM_IMAGE_ACQUISITION_MODE	TSI_TYPE_UN32: Obsolete, do not use
TSI_PARAM_NAMED_VALUE	TSI_TYPE_TEXT: For internal use only.
TSI_PARAM_TAPS_INDEX	TSI_TYPE_UN32: Used to select or query the current tap setting.
TSI_PARAM_TAPS_VALUE	TSI_TYPE_UN32: Used to query the number of taps for the current tap setting (See TSI_PARAM_TAPS_INDEX above).
TSI_PARAM_RESERVED_1	Reserved for internal use.
TSI_PARAM_RESERVED_2	Reserved for internal use.
TSI_PARAM_RESERVED_3	Reserved for internal use.
TSI_PARAM_RESERVED_4	Reserved for internal use.
TSI_PARAM_GLOBAL_CAMERA_NAME	TSI_TYPE_TEXT: Returns the camera name stored in the camera thats visible to users on the network and software loading the camera.
TSI_PARAM_CDS_GAIN_VALUE	TSI_TYPE_UN32: For expert use only, contact us for details.
TSI_PARAM_PIXEL_SIZE	TSI_TYPE_FP: Size of the pixel width in microns (pixel assumed square).
TSI_PARAM_BITS_PER_PIXEL	TSI_TYPE_UN32: Number of bits per pixel in raw data.
TSI_PARAM_BYTES_PER_PIXEL	TSI_TYPE_UN32: Number of bytes per pixel (1 or 2).
TSI_PARAM_READOUT_TIME	TSI_TYPE_FP: Readout time of the camera in seconds.
TSI_PARAM_HW_TRIGGER_ACTIVE	TSI_TYPE_UN32: Gets or sets the hardware active trigger state. 1 = Active, 0 = Not Active.
TSI_PARAM_HW_TRIG_SOURCE	TSI_TYPE_UN32 (enum): Gets or sets the HW trigger source.
TSI_PARAM_HW_TRIG_POLARITY	TSI_TYPE_UN32 (enum): Gets or sets the HW trigger polarity.
TSI_PARAM_TAP_BALANCE_ENABLE	TSI_TYPE_UN32 : Gets or sets the tap balance enable. 1 = Active, 0 = Not Active.

2.5. TsiImage Class

2.5.1. DATA MEMBERS

```

unsigned int m_Width;           // Width of image in pixels.
unsigned int m_Height;          // Height of image in pixels.
unsigned int m_BitsPerPixel;    // The number of significant
                                // bits per pixel in the
                                // pixel data.

unsigned int m_BytesPerPixel;    // The number of bytes
                                // consumed by a pixel.
unsigned int m_SizeInPixels;     // Size of image in pixels.
unsigned int m_SizeInBytes;      // Size of image in bytes.
unsigned int m_XBin;             // Horizontal binning value.
unsigned int m_VBin;             // Vertical binning value.
unsigned int m_ROI[4];           // The region of interest
                                // (subimage) that the pixels
                                // were gathered from.
                                // Format: [x0, y0, x1, y1]
unsigned int m_ExposureTime_ms;  // Exposure time in
                                // milliseconds.
unsigned int m_FrameNumber;      // Frame number returned from
                                // frame grabber.

union
{
    void          *vptr;
    char          *i8;
    unsigned char  *ui8;
    short         *i16;
    unsigned short *ui16;
    unsigned int   *ui32;
} m_PixelData;

```

2.5.2. METHODS

TsiImage::Copy

```
bool Copy(TsiImage *src)
```

Copies an existing TsiImage's data.

TsiImage::Clone

```
TsiImage* Clone(TsiImage *src)
```

Creates a new duplicate of the specified TsiImage.

2.6. TsiUtil Class**2.6.1. METHODS*****TsiUtil:WriteImageToTIFF***

```
bool WriteImageToTIFF (char *file_name, TsiImage *image)
```

Saves the TsiImage with the specified name into a 16-bit Tiff image.

TsiUtil:WriteImageToPNG

```
bool WriteImageToPNG (char *file_name, TsiImage *image)
```

Saves the TsiImage with the specified name into a 16-bit PNG image.

Chapter 3 Thorlabs Worldwide Contacts

USA, Canada, and South America

Thorlabs, Inc.
56 Sparta Avenue
Newton, NJ 07860
USA
Tel: 973-579-7227
Fax: 973-300-3600
www.thorlabs.com
www.thorlabs.us (West Coast)
Email: sales@thorlabs.com
Support: techsupport@thorlabs.com

Europe

Thorlabs GmbH
Hans-Böckler-Str. 6
85221 Dachau
Germany
Tel: +49-(0)8131-5956-0
Fax: +49-(0)8131-5956-99
www.thorlabs.de
Email: europe@thorlabs.com

France

Thorlabs SAS
109, rue des Côtes
78600 Maisons-Laffitte
France
Tel: +33 (0) 970 444 844
Fax: +33 (0) 825 744 800
www.thorlabs.com
Email: sales.fr@thorlabs.com

Japan

Thorlabs Japan, Inc.
Higashi-Ikebukuro Q Building, 1F
2-23-2, Higashi-Ikebukuro,
Toshima-ku, Tokyo 170-0013
Japan
Tel: +81-3-5979-8889
Fax: +81-3-5979-7285
www.thorlabs.jp
Email: sales@thorlabs.jp

UK and Ireland

Thorlabs Ltd.
1 Saint Thomas Place, Ely
Cambridgeshire CB7 4EX
Great Britain
Tel: +44 (0)1353-654440
Fax: +44 (0)1353-654444
www.thorlabs.com
Email: sales.uk@thorlabs.com
Support: techsupport.uk@thorlabs.com

Scandinavia

Thorlabs Sweden AB
Mölnadalsvägen 3
412 63 Göteborg
Sweden
Tel: +46-31-733-30-00
Fax: +46-31-703-40-45
www.thorlabs.com
Email: scandinavia@thorlabs.com

China

Thorlabs China
Room A101, No. 100
Lane 2891, South Qilianshan Road
Putuo District
Shanghai
China
Tel: +86 (0)21-60561122
Fax: +86 (0)21-32513480
www.thorlabs.hk
Email: chinasales@thorlabs.com



THORLABS
www.thorlabs.com

